

Dilution Refrigerator with cryogenic cycle circulation and adsorption pumping

Vladimir Mikheev, Oxford

***Compact self contained Cryofree continuous Dilution refrigerator down to 30mK
(Oxford Instruments 2006-2008)***

List of participants:


**Alvin Adams ,
Tim Foster
Rod Bateman,
Paul Noonan
Vladimir Mikheev**



50 Years DR

Size Evolution for DR at OI


70s'



80s'




2000



50 Years DR

The way to further miniaturization of Dilution refrigerators is in Developing Cryogenic cycle of He3 circulation. Why?



50 Years DR


Dilution Refrigerators with cryogenic circulation of He3 at Oxford Instruments since 1997

- The first example is **AST Minisorb** top loaded in liquid He4.
- Compact DR, --16 mK, quick turn around (<3 hrs from RT)

AST Minisorb Issues:

- Two cryo-valves ,
- A number of high value impedances,
- 1K Pot external pump.

- No external pumps,*
- No cryo-valves,*
- No impedances*
- New technology - Cryofree Technology*




50 Years DR


Cryogenic cycle of He3 circulation

How to arrange
Cryogenic circulation?

Heat Pipe effect



Powerful continuous He3 fridge



Continuous He3-He4 Dilution refrigerator



50 Years DR

Pre-History:

- **H.London, G. Clarke, E.Mendoza Phys.Rev. 128,1992, 1962**
- **V Edelman , Cryogenics 12.385.1972 (T~60mK @ Single shot operation)**
- **A. Aroyan, R.Amamchan et al (preprint, Erevan) T~0.2K, attempt of continuous**
- **V Sobolev, V Syvokon et al.1992, Cryogenics ICEC Supplement (12mK continuous)**
- **Yu Bunkov, J Niekky et al.. ~1992 QFS (95-96)(12mK, continuous)**
- **V Edelman (with Institute Photonic, Berlin) 2006, (T~50mK, Single shot)**



50 Years DR

The two main functions at the operation DR refrigerators:

1. Condensation of gas
2. Circulation of gas or liquid.


The key thing is to produce a liquid He4 or He3, (liquefaction)

What is the best way of He3 liquefaction? 1Kpot?




50 Years DR

Condensation on the cold surface (1K Pot method)



50 Years DR

- **Condensation on the cold surface requires:**
- **(1K Pot + $N/V_{or}(Z)$ + External Pump)**



50 Years DR

- Liquefaction of He3 by adiabatic expansion gas is easier
- Creating the “rain” of liquid He3 or He4 at the pressure drop



50 Years DR



**Adiabatic expansion in atmosphere
Creating the “rain”**

50 Years DR



50 Years DR

$$PV = mRT \quad \text{Ideal Gas}$$

Reversible adiabatic process

$$P_1 V_1 T_1 \xrightarrow{S = \text{const}} P_2 V_2 T_2$$

$$\Delta S = c_v \ln\left(\frac{T_2}{T_1}\right) + R \ln\left(\frac{V_2}{V_1}\right) = 0$$

$$\left(\frac{T_2}{T_1}\right) \left(\frac{V_2}{V_1}\right)^{R/c_v} = 1$$


$$\left(\frac{T_2}{T_1}\right) = \left(\frac{P_2}{P_1}\right)^{(\gamma-1)/\gamma} \quad \gamma = c_p/c_v = 5/3$$

$$T_2 = T_1 \left(\frac{P_2}{P_1}\right)^{0.4}$$

Example: $T_1 = 4.2\text{K}$; $P_2/P_1 = 1/3$

$$T_2 = 2.7\text{K} < T_{\text{condens}}$$


He3




50 Years DR

From mechanical expanders to non-mechanical way of expansion


Piston expander



Expansion by cooling of reservoir





Single expansion by F.Simon, 1932
(method of L.Cailletet, 1877)
Troitskii (1980)



50 Years DR

- Cryomech PT 403 Cold Head
- Electrical Connectors
- ^3He Cold Dumps
- PTR first stage thermal link & shield connection
- Adsorption pump
- Heat switch
- PTR second stage thermal link & shield connection
- He^3 Pot and sample region at $<280\text{mK}$ for 2-3 days



© Oxford Instruments

General view of He^3 system with PT403,
1 phase, air cooled compressor



50 Years DR

Novel self contained continuously operating Dilution Refrigerator down to mK range (Ox.Inst. 2006)

Key features:



- **Continuous Cryogenic Circulation – CCDR**
- **No External Pumps**
- **No Gas Handling System, - whole gas mixture is inside the Fridge**
- **No liquid cryogen is used, even no cold traps and cold spots on the system**
- **Low Power Single Phase Air cooled compressor (PTR<250mW@4K)**
- **Push button operation, with computer control.**




50 Years DR

Principal Diagram of CCDR

Continuous operation




50 Years DR




50 Years DR

Cryofree Continuous operation of He3 0.4K Stage for Triton-DR



50 Years DR

6 days of continuous operation
Cryogen free Dilution refrigerator Triton™DR



50 Years DR

Specification

CCDR specification:

Base Temperature 20-50mK


Cooling power @ 350 mK

for thermal dumping all heat leaks.....> 200 μ W

50 Years DR

Cryogenic Cycle of He3 Circulation presents new type of self contained Dilution refrigerators

- Cryogen Free (no liquid cryogen service, no cold traps)
- Continuous operation, at temperature down to 20mK so far
- Very compact, no external pumps, no gas handling system
- Simple T-control, Femto power MC temperature control
- True turnkey; plug in the wall, stand-alone control system with PC user interface via Ethernet (done from Oxford to Goteborg University)
- Low vibration (due to small PTR)
- Large sample space and easy sample change, creates opportunity to interface cryostat to a wide range of applications and opens ULT to a wider audience in the scientific community
- ULT will be re-invented and de-mystified...



50 Years DR



Thank you for your
attention

Any questions?

50 Years of Dilution Refrigeration!

50 Years DR

Demand of Cooling Power of DR

2 options

```
graph TD; A[2 options] --> B[Standard solution]; A --> C[CCDR solution]; B --> D[High cooling power DR at Mixing Chamber]; C --> E[Additional powerful ULT below platform 0.8K and reduced demand on cooling power at Mixing Chamber];
```



Standard solution

CCDR solution

**High cooling power DR at
Mixing Chamber**

**Additional powerful ULT below
platform 0.8K and reduced
demand on cooling power at
Mixing Chamber**

50 Years DR



$$\frac{\dot{Q}_{\text{standard}}}{\dot{Q}_{\text{CCDR}}} \approx \left[\frac{0.8K}{0.35K} \right]^{\geq 2}$$

> 5 for metals

>25 for Insulators

50 Years DR

Cooling Power of DR

Required for:

- MC temperature low T, stability, heat withdrawal from customer's experimental interventions, large sample
- Dumping all heat leaks coming from services : wiring, cables, Drivers, optical access lines, additional instrumentation (SQUID)

Base Temperature is Lowest temperature on MC at

$$Q_{\text{heat load}} = Q_{\text{cooling power on MC}}$$

50 Years DR

Heat Load

Heat transfer through materials at Low temperatures

$$\dot{Q} = \frac{S}{L} \int_{T_{MC}}^{T_d} k(T) dT$$

$$T_{MC} \ll T_d$$

Metals, $k \sim T$
Conduction electrons,

$$\dot{Q} = \frac{S}{2L} T_2^2 dT$$

Insulators/Phonons
 $K \sim T^3$; @ $T < \Theta_D/10$

$$\dot{Q} = \frac{Sb}{4L} T_2^4 dT$$

50 Years DR

Control System

Triton™-DR

- **FPGA (field programmable gate array) based algorithms**
- **Automated cycle, stand-alone operation**
- **Femtopower™ technology for MC T-sensor**
- **User interface via ethernet from anywhere**
- **Password security option**